


## Dual Precision Syringe Pump

Renovo made an easy to operate precise and reliable syringe pump with feature options for your lab requirements. You can use the device only with one hand for lab practicality.

Easy to use precise calibration settings designed to ensure you have full control of the device from the moment you start to use.

Optional precision heater control has been added for your special needs arising in laboratories especially on bio-materials and cell culture operations.

Easy slide and easy syringe locking equipment is designed for you to use it practically even when one hand is full.



### Heated Syringe

Heated Syringe Bed allows up to 60 C heated injections for project requirements.



### Multiple Pumps

Up to 2 optional syringe pumps can be operated sequentially stop pause reset and work together.



### Precise Control

Precision options down to 0.1µl per hour and µl based volume setting allows you to precise control the injection



### Easy to Use

Standalone button panel and Easy Slide pump pusher makes easier start stop refill operationsn



### One Hand Use

Sliders and syringe locks can be used with only one hand making device more useful

### Fast Calibration

Precise fast calibration allows user to set up process faster than any other device

## TECHNICAL SPECIFICATIONS

<b>Amount of Pumps</b>	1 or 2	<b>Step Resolution</b>	1mm/1200 step
<b>Suitable Syringes</b>	0.5µl to 150 ml	<b>Force</b>	100N
<b>Heater</b>	Optional Up to 60°C	<b>Synchronised Start</b>	2 Syringes
<b>Suitable Power</b>	220v	<b>Syringe Pusher</b>	One hand easy slide
<b>End Stop</b>	Yes	<b>Syringe Lock</b>	One hand easy lock
<b>Max Speed</b>	200mm/min	<b>Calibrate Mode</b>	Precise control
<b>Dimensions</b>	420x330x250mm	<b>Weight</b>	11.2kg